

## GC sample preparation protocol

For GC-TOF-MS sample extracts (BR: 200  $\mu\text{L}$ ; OB: 150  $\mu\text{L}$ ) containing internal standards were evaporated under a stream of nitrogen gas at 37°C. Derivatisation was achieved via oximation and silylation. To oximate the samples, methoxyamine hydrochloride (20 mg/mL in pyridine, BR: 50  $\mu\text{L}$ ; OB: 25  $\mu\text{L}$ ) was added and the samples were vortex mixed ( $\sim 1$  min) to dissolve the dried compounds. The samples were then incubated at 60 °C for 1 hour. Thereafter, the samples were silylated by adding BSTFA (BR: 50  $\mu\text{L}$ ; OB: 25  $\mu\text{L}$ ) containing 1 % (v/v) TMCS and incubated at 60 °C for 1 hour. For the BR samples, the final volume (100  $\mu\text{L}$ ) was then transferred to 250  $\mu\text{L}$  pulled point glass inserts. Due to the small volumes used during derivatisation, OB samples (50  $\mu\text{L}$ ) were derivatised in microvials designed to maximise recovery for small sample volumes. Finally, each vial was loaded onto an Agilent© 7693 auto sampler for GC-TOF-MS analysis.