

<b>West Coast Metabolomics Center</b>	<b>SOP Standard Operating Procedure</b>	<b>page 1 of 3</b>
<b>date:</b> 06/24/2013	<b>Sample preparation of blood plasma or serum samples for lipidomic analysis</b>	Code no.: PlasmaLipidExtraction 06242013

Issued: 06/24/2013 Valid from: 06/24/2013	Validity area: UC Davis Genome Center, Metabolomics Core and Research Laboratories
Responsible: Bill Wikoff	Secondary: Carlos Leon, Brian DeFelice
This SOP supersedes: Sample preparation of blood plasma or serum samples for lipidomics analysis 12/17/2012	approved: Oliver Fiehn

### **Sample preparation of blood plasma/serum samples for lipidomic analysis**

#### **1. Purpose**

This SOP describes sample extraction and sample preparation for lipid profiling by liquid chromatography/quadrupole time-of-flight mass spectrometry (LC-QTOF) or nanoelectrospray ion trap-FTICR MS.

#### **2. References**

Matyash V, Liebisch G, Kurzchalia TV, Shevchenko A and Schwudke D (2008) Lipid extraction by methyl-*tert*-butyl ether for high-throughput lipidomics. *J Lip Res* 2008, 49: 1137-1146

#### **3. Starting material**

Plasma/serum: 20 µl sample volume or aliquot (see Aliquoting TEDDY Plasma SOP)

Control: Pooled Sterile Human Plasma in Sodium Citrate from Rockland Institute (aka "Citrate Plasma"). Citrate Plasma 1 mL aliquots are stored in the Revco -80°C freezer at 2-6.

#### **4. Equipment**

- Centrifuge (Eppendorf 5415 D)
- Calibrated pipettes 20–200 µL and 100-1000µl
- Eppendorf tubes 1.5 mL, uncolored (Cat. No. 022363204)
- ThermoElectron Neslab RTE 740 cooling bath at -20°C
- MiniVortexer (VWR) 58816-121
- Orbital Mixing Chilling/Heating Plate (Torrey Pines Scientific Instruments)
- Speed vacuum concentration system (Labconco Centrivap cold trap)
- Eppendorf tips for organic solvents such as acetonitrile, methanol, and MTBE
- Glass Amber Vials: National Scientific (C4000-2W)
- Glass Inserts: Supelco 27400-U
- Blue Tops for Vials: Agilent (5182-0717)

#### **5. Chemicals:**

- Crushed ice
- Nitrogen line with pipette tip
- Pure LC/MS Grade Water (Fisher Optima W6-4)
- MTBE: Sigma, Chromasolv 99.8% for HPLC 100mL (smallest available) (34875-100mL)
- Methanol: J.T. Baker LC/MS Grade (9830-03)
- CUDA (12-[[[(cyclohexylamino)carbonyl]amino]-dodecanoic acid) from Cayman Chemical Item Number 10007923
- Ethanol absolute (Sigma-Aldrich)

#### **6. Sample Preparation**

Extraction solvents

<b>West Coast Metabolomics Center</b>	<b>SOP Standard Operating Procedure</b>	<b>page 2 of 3</b>
<b>date:</b> 06/24/2013	<b>Sample preparation of blood plasma or serum samples for lipidomic analysis</b>	Code no.: PlasmaLipidExtraction 06242013

- Degas both “MeOH with QCmix” and MTBE with Cholesterol Ester 22:1 (see SOP “QC mix for LC-MS lipid analysis”) by sonication
- Store solvents in the –20°C freezer to pre-chill

#### Homogenization and extraction

- Thaw each 20 µL plasma aliquot at room temperature (see Aliquoting TEDDY samples SOP). Once thawed (~10min) place liquid plasma samples on ice.
- Add 225 µL cold “MeOH with QC mix” (see SOP “QC mix for LC-MS lipid analysis”). Keep MeOH on ice during extraction
- Vortex each sample for 10s, keeping the rest on ice during all the extraction.
- Add 750 µL of cold MTBE with 22:1 CE, keep MTBE on ice during extraction
- Vortex for 10s
- Shake for 6min at 4°C in the orbital mixer.
- Add 188 µL room temperature LC/MS grade water.
- Vortex for 20 s
- Centrifuge for 2 min @ 14,000 rcf (12300 rpm)
- Remove supernatant, splitting into two aliquots of 350 µL, keeping one at –20°C for backup
- Dry samples to complete dryness in the speed vacuum concentration system

#### Preparation of resuspension solvent

- Weight 2mg of CUDA and dissolve in 2mL of ethanol. This is the CUDA stock solution (1mg/mL).
- Prepare 120 mL of MeOH:Toluene 90:10 (108 mL MeOH+12 mL Toluene). Degas by sonication for 5min
- Sonicate CUDA stock solution (1mg/mL) for 5 min, then dilute by adding 50µL (1mg/mL CUDA) in 950 µL MeOH:Toluene 90:10. This is CUDA 0.05mg/mL.
- Dilute 100 µL CUDA 0.05mg/mL in 9.9 mL MeOH:Toluene 90:10. This is CUDA 0.5 µg/mL.
- Dilute 10mL of CUDA 0.5 µg/mL in 90 mL MeOH;Toluene 90;10. This is the final resuspension solvent, MeOH:Toluene 90:10 with CUDA 50ng/mL.
- Store the solutions in the freezer at -20°C until use.

#### Preparation for analysis

- Re-suspend dry samples in 108.6 µL MeOH:Toluene 90:10 with CUDA (50ng/mL), 24 samples at a time, degassed using the above method.
- After adding MeOH:Toluene to all 24 samples, vortex the rack on low speed as to not “whip” air into the liquid for 20 seconds.
- Then Sonicate at RT for 5 min, and centrifuge samples for 2 min at max speed.
- Transfer 50 µL to two separate amber glass vial with micro-insert. Cap vials with Agilent blue top.
- Use independent vials for positive and negative mode acquisitions. Use the QTOF 6530 for positive ion analysis and the QTOF 6550 for negative analysis.

#### **7. Problems**

To prevent contamination disposable material is used. To prevent inhalation of toxic ether vapor, use fume hood during lipid extraction.

<b>West Coast Metabolomics Center</b>	<b>SOP Standard Operating Procedure</b>	<b>page 3 of 3</b>
<b>date:</b> 06/24/2013	<b>Sample preparation of blood plasma or serum samples for lipidomic analysis</b>	Code no.: PlasmaLipidExtraction 06242013

#### 8. Quality assurance

- For each sequence of sample extractions, perform one blank negative control extraction by applying the total procedure (i.e. all materials and plastic ware) without biological sample.
- Use TEDDY citrate plasma standard aliquot per 10 authentic subject samples as control.
- Prepare at least six NIST plasma extracts in the same manner as positive controls

#### 9. Disposal of waste

- Collect all chemicals in appropriate bottles and follow the disposal rules.
- Collect residual plasma / serum samples in specifically designed red “biohazard” waste bags.

### Final Protocol

