

# HepG2\_FN3K Sample Preparation Protocol

## Description

Protocol adapted from SECIM SOP: NMR sample prep of cells, effective 11/24/2014

Original protocol cited for approx 6 million cells.

## Materials and Reagents

PBS buffer (15 mL x # of plates)

Aspirator

Dry ice

Liquid nitrogen

Cell scraper (one per plate)

Refrigerated centrifuge

Extraction solvent (Methanol:Water, 80:20 (v/v) UPLC or LC-MS grade)

Speed-vac concentrator

3mm NMR tubes

100 mM sodium phosphate buffer in D<sub>2</sub>O, pH 7.4

Includes: 1/3 mM Sodium trimethylsilylpropanesulfonate (DSS-D6 )

3 mM sodium azide

## Sample Collection and Extraction Procedure

Rapid quenching of metabolic reactions is required to minimize the occurrence of unwanted metabolic stress and 'run-on' reactions.

### Steps

1	Wash the monolayer three times with 5 ml <u>cold</u> PBS.
2	Add <u>ice cold</u> extraction solvent to plate, scrape cells, and transfer to a tube. (10 cm plate: 1.0 ml extraction solvent)
3	Freeze tubes immediately in liquid nitrogen and store at -80°C until extraction.
4	Thaw tubes on ice. (Approx 1 hr)
5	Vortex tubes well twice for one minute intervals.
6	Pellet the debris by centrifugation on highest speed for 15 min at 4°C.
7	Transfer the supernatant to a clean tube and label.
8	Keep debris tube as well for protein concentration determination/sequential extraction.
9	Evaporate the solvent from extracts using speed-vac instrument until dry, approx. 4-10 hrs (Labconco).
10	Store dried samples at -80°C until NMR/mass spec analysis.

Steps 1-3 Performed by Sami and Kannan Lab.

Harvested cells transferred on dry ice to Edison Lab 2/8/18.

Rest of extraction procedure performed by MBC in Edison Lab.

---

## NMR Preparation

Dried extracts resuspended in 200  $\mu$ L sodium phosphate NMR buffer in randomized run order. Buffer blanks loaded directly into tubes before first sample and after last sample resuspended. All steps performed on ice/4°C.

1. After addition of buffer, samples vortexed for 1 min
2. Spun down in microcentrifuge for 15 sec.
3. 180  $\mu$ L transferred into 3mm NMR tubes.
4. 15  $\mu$ L of remaining sample combined to form internal pooled samples, loaded last.
5. NMR tubes kept on ice/4°C until transfer to SampleJet for analysis