

Extraction Protocol: Liver

Platforms: CSH-QTOF (Lipidomics)/ GCTOF (Primary Metabolites)/ Proteomics/ Oxylipins

Amount of Tissue Provided: ~100-150mg

- 1) Keep Specimen on Dry Ice
- 2) Transfer Tissue Contents into a new 1.5mL labeled eppendorf tube; keep on dry ice at all times
- 3) Add three (3) 3mm metal grinding balls to each sample; store in -80C for 10minutes
- 4) Homogenize the entire tissue to fine powder using genogrinder; make sure that the metal grinding balls are ice-cold prior to homogenization (step 3)
- 5) Upon completion of homogenization, keep samples on dry ice
- 6) Weight out **two (2) aliquots**: a ~**5mg** aliquot for *CSH_lipidomics* and a ~**4mg** aliquot for *Primary Metabolites by GCTOF*
 - a. Record the exact weight weighed out for each sample
 - b. Keep all samples on dry ice**
- 7) KEEP remaining tissue specimen (>90mg) for analysis of Oxylipins (store in -80C)
- 8) Analysis of Primary Metabolites (GCTOFMS)**
 - a. Add 1mL of ice-cold “degassed” 3:3:2 ACN/IPA/H₂O
 - b. Vortex for 10seconds
 - c. Shake on shaker for 20min at -4C
 - d. Centrifuge the samples for 2min at 14,000 rcf
 - e. Transfer two (2) 500µL aliquots to new 1.5mL eppendorf tubes; one for backup the other to be dried to dryness using the SpeedVac
 - f. **IMPORTANT:** The precipitated protein will be used for analysis of the proteome, DO NOT DISCARD THESE; Place these in a separate labeled box and store in -20C
 - g. Keep all samples on ice during extraction period
 - h. Dry down one (1) 500µL aliquot to complete dryness
 - i. Perform cleanup on dried aliquot using 500µL of 50/50 v/v ACN/H₂O
 - j. Transfer supernatant and dry to completeness
 - k. Submit for Derivatization
- 9) Analysis of Complex Lipids (LCQTOF)**
 - a. Add 225µL of ice-cold “degassed” MeOH containing “ISTD mixture” to homogenized 5mg aliquot
 - b. Vortex for 10 seconds
 - c. Add 750µL of ice-cold “degassed” MTBE containing 22:1 CE ISTD
 - d. Vortex for 10 seconds
 - e. Shake on Orbital Mixer for 6min at 4C
 - f. Add 188µL of room temperature H₂O
 - g. Vortex for 20 seconds
 - h. Centrifuge for 2min at 14,000 rcf
 - i. Transfer two (2) aliquots of 350µL of top layer, one for backup stored in -20C, the other for analysis
 - j. Keep bottom layer and store in -20C
 - k. Dry down one (1) 350µL aliquot to dryness using the Speedvac
 - l. Resuspend samples in 108.6µL of 50ng/mL CUDA
 - m. Vortex and sonicate for 5minutes
 - n. Centrifuge for 2min at 14,000 rcf
 - o. Transfer 90µL to an amber vial with micro-insert (non-diluted)

- p. Transfer 10 μ L to a new 1.5mL eppendorf tube, dilute 20X with 50ng/mL CUDA in 90:10 MeOH:Toluene (10 μ L + 190 μ L CUDA) and transfer 100 μ L to amber vial with micro-insert (diluted for TGs)
 - i. The dilution is based off previous experiences with liver samples (as described by Ingrid)