Sample preparation for sulfonate measurements

Solvents/Chemicals:

Acetonitrile was purchased from Sigma-Aldrich (St. Louis, MO, USA). All solvents for mass spectrometry were of analytical grade purity. Water (resistivity 18.2 MΩcm) was purified using a Milli-Q-System (Millipore, Milford, MA, USA). Formic acid was purchased from Honeywell Fluka (Muskrgon, MI, USA). Cysteic acid and isethionic acid were purchased from Sigma Aldrich (St. Louis, MO, USA). 2,3-dihydroxypropane-1-sulfonate and sulfoquinovose were purchased from MCAT GmbH (Donaueschingen, Germany). Taurine was purchased from Roth GmbH(Karlruhe, Germany). Ti(III) nitrilotriacetate (Sigma-Aldrich, Germany)

Anoxic phosphate-buffered saline (PBS): 8.5~g / I NaCl, 0.3~g / I KH2PO4, Na2HPO4 0.6~g / I, 0.1~g / I Bacteriological Peptone, 1~mg / I Resazurin, 40~mM Sodium DL-lactate , 40~mM Sodium formate pH 7.0, N2/CO2 (80/20, v / v) as gas phase, autoclaved at 121° C for 15~min. all compounds from Fluka ,Muskrgon, MI, USA.

DS-medium: 19 mM NH4Cl, 17 mM NaCl, 2 mM MgCl2 x 6 H2O, 7 mM KCl, 0.3 mM CaCl2 x 2 H2O, 1 mM K2HPO4, 40 mM sodium DL-lactate, 40 mM sodium formate, 3.5 mg/l yeast extract, 1 ml/l selenite-tungstate solution 1 ml/l trace element solution (see Table 5), 1.2 μ M 1,4-naphthoquinone and 2 μ M resazurin. The medium was adjusted to pH 7.4, gas flashed with N2/CO2. all compounds from Fluka ,Muskrgon, MI, USA.

Solution	Components with concentration
Trace-element solution	10 ml/l HCl
	1.5 g/l FeCl ₂ x 4 H ₂ O
	70 mg/l ZnCl ₂
	100 mg/l MnCl ₂ x 4 H ₂ O
	6 mg/l H₃BO₃
	190 mg/l CoCl ₂ x 6 H ₂ O
	2 mg/l CuCl ₂ x 2 H ₂ O
	24 mg/l NiCl ₂ x 6 H ₂ O
	36 mg/l Na ₂ MoO ₄ x 2 H ₂ O
Selenite-tungstate solution	500 mg/l NaOH
	3 mg/l Na ₂ SeO ₃ x 5 H ₂ O
	4 mg/l Na ₂ WO ₄ x 2 H ₂ O

Seven-vitamin solution	100 mg/l vitamin B ₁₂
	80 mg/l p-amino benzoic acid
	20 mg/l D (+)-biotin
	200 mg/l nicotinic acid
	100 mg/l calcium pantothenate
	300 mg/l pyridoxine hydrochloride
	200 mg/l thiamine hydrochloride x 2 H ₂ O
Ti(III) nitrilotriacetate solution	19.2 g/l nitrilotriacetic acid diluted in anoxic
	distilled water, pH of 9 adjusted with NaOH
	19.2 ml 20% TiCl₃ (Acros)
	pH of 7adjusted with Na₂CO₃ (80 g/l)

Work steps calibration samples

One g of the human fecal sample was homogenized in anoxic phosphate-buffered saline (PBS, pH 7.0, Table 4) by vortexing with glass beads (c. 3 mm; Roth, Germany) to yield a 10% fecal suspension. The fecal suspension was further diluted to 1% in a Hungate tube containing anoxic PBS supplemented with 3.18 mM of sterile filtered Ti(III) nitrilotriacetate (Sigma-Aldrich, Germany) as reductant. Subsequently, fecal slurry aliquotes were centrifuged (14.000 x g, 4° C, 5 min) and the supernatant frozen until further processing. Fecal slurries were , fecal supernatants were serially diluted using 50% acetonitrile and 0.1% formic acid in water to obtain final dilution of 1:10,000 and spiked with cysteate, sulfoquinovose (SQ), 2,3-dihydroxypropane-1-sulfonate (DHPS), taurine and isethionate for calibration. 500 μ l of each sample were placed in glass vials (Wicom, Heppenheim, Germany) and stored at - 80 °C until LC-MS/MS analysis.

Work steps for application pilot study

For each 10 ml incubation, DS medium supplemented with either taurine (20 mM) or sulfoquinovose (SQ, 4 mM) was inoculated with the SIHUMI bacteria or the SIHUMI consortium and B. wadsworthia. Control incubations of medium containing only sulfonates or bacteria were included. The incubations were conducted under anoxic conditions at 37 °C in duplicates. Samples (600 μ l) were withdrawn after 0, 3, 24 and 48 h for quantification of sulfoquinovose (SQ), 2,3-dihydroxypropane-1-sulfonate (DHPS), taurine and isethionate. For the analysis of sulfonates by LC-MS/MS-MRM, aliquots of 250 μ l for each sample were centrifuged (14 000 x g, 4 °C, 5 min) and 50 μ l of the supernatant stored at – 20 °C until further processing. Samples were thawed, centrifuged at 18 000 x g at RT for 2 min and the supernatant was diluted 1:10 000 in 50% aqueous acetonitrile. Subsequently, 500 μ l of each sample were placed in glass vials (Wicom, Heppenheim, Germany) and stored at - 80 °C until LC-MS/MS analysis.