

Precellys Homogenization and Quenching Procedure

For every tissue sample we use:

2 ml 100 % acetonitrile

1.3 ml dH₂O

0.2 ml 0.2 mM Tris-Cl pH8.0 (40 nmoles of Tris)

Total volume = 3.5 ml

To avoid mixing up 60%AcN with Tris and 60% AcN without Tris, my suggestion is 1

We can make master mix for 20 samples as follows:

Master Mix:

2 ml x 20 = 40 ml 100% acetonitrile

1.3 ml x 20 = 26 ml dH₂O

0.2 ml x 20 = 4 ml 0.2 mM Tris-Cl pH 8.0 (40 nmoles Tris)

Total V = 70 ml

Slice Pre-treatment:

- * Add 100 ml master mix to snap-cap tube (depending on the slice size and shape)
- * Use microscissors to cut the slices inside the snap-cap tube
- * Use pestle to smash the pieces against the tube wall (BE VERY CAREFUL)
- * Use microscissors again to cut the slices inside the snap-cap tube
- * Use pestle again to smash the pieces against the tube wall (BE VERY CAREFUL)
- * Wash the pestle with 100 ml master mix
- * Wash the scissors with 100 ml master mix
- * Transfer sample to Precellys tube containing beads using 1 ml pipet tips cut at the tip
- * Wash snap-cap tube and pipet tip with 2 x 100 ml master mix. KEEP the tip.

In this way, about 500 ml of pre-cut slices will be further homogenized using Precellys:

Precellys:

Add ice to coil and wait till air T goes down to 16 °C.

Use 3 cycles of 5000 rpm - 5 seconds interval in between runs.

At this point, you have 3 ml of master mix to use for transferring the homogenized

Rinse the beads and tube with 3 x 1.0 ml of master mix, transferring to the 15 ml tube

to prepare a master mix with all 3 components described above.

you may need to use another 50 or 100 ml - if this is the case, remember that the f

the tip end with razor blade to facilitate transfer to small tissue chunks.

allys.

d/quenched sample to 15 ml conical tube and washing the tube with beads and the
tube.

final volume in the Precellys tube with beads should be around 500 ul)

⇒ pipet tip you saved before.