

20170322 2HG extraction procedures for cell matrix

Samples

Castro-Lowenstein, EX00705, N=24 cells, LOC 0000
, EX00 , N= cells, LOC0000

Extraction Solvent

_____ mL of 7:2:1 (Methanol: H₂O: Chloroform) + _____ μ L of 100 μ M ¹³C Hydroxyglutaric acid [Final 2 μ M]

DATAN: diacetyl-L-tartaric anhydride- 50mg/mL in Dichloromethane-**acetic acid** (4:1- 5 parts)

Standard Mix stock (STD) and Internal Standard stock (IS)

Hydroxyglutaric acid ¹³C 100 μ M

Standard Preparation (STD Mix)--LCMS

STD [Final- μ M]	MeOH (μ L)	ChCl ₃ (μ L)	H ₂ O (μ L)	2HG STD mix (μ L) of 25 μ M D + L	Hydroxyglutaric Acid 100 μ M
STD 0 [0 μ M]	140	20	40	0	4
STD 1 [0.06 μ M]	140	20	39.5	.5	4
STD 2 [0.125 μ M]	140	20	39	1	4
STD 3 [0.625 μ M]	140	20	35	5	4
STD 4 [1.25 μ M]	140	20	30	10	4
STD 5 [2.5 μ M]	140	20	20	20	4
STD 6 [5 μ M]	140	20	0	40	4
STD A	100	20	0	80 D-2HG	4
STD B	100	20	0	80 L-2HG	4

Sample Preparation

- 1) Place all samples on wet ice until extraction procedure has been completed.
- 2) While over wet ice **add 0.6ml Extraction solvent** containing Internal Standards (ISs) samples.
- 3) Sonicate with probe for 5-10secs set to 20% Duty cycle and 2 output, repeat vortex.
- 4) Keep eppendorf tubes at **4°C for 10min** to allow complete extraction, remove from 4°C repeat vortex.
- 5) **Centrifuge** all tubes at 14,000RPM for 10min in 4°C.
- 6) **Create a pooled** sample by transferring equal volumes from each sample to an autosampler vial.
- 7) **Transfer 200 μ L** of supernatant to an autosampler vial take to dryness using speed vac at 55° Aqueous for ~1 hr.
- 8) **Add 50 μ L** of 50mg/mL DATAN, cap and **incubate** at **75°C for 30 min**.
- 9) **Cool vials and dry** by continuous N₂ flow at RT°C for ~ 1 hr.
- 10) **Reconstitute** samples in 100 μ L of LC grade H₂O, vortex, and transfer to insert.

LC-MS Analysis

QQQ- method: 2HG 6490 in CR methods folder

Column: Waters HSS T3 C18 50mm column.

MP A: 2mM Ammonium Formate in H₂O, pH ~3.3-4.01, adjust with LC-MS grade Formic Acid

MP B: 100% ACN with 0.1% Formic Acid

$D = 363.2 \rightarrow 147.2$

Notes/Observations